

3D-Printed Implants: How Additive Manufacturing Is Personalizing Orthopedic Medicine

Patient-specific implants printed from titanium and PEEK are transforming bone repair — here is the materials science making it possible.

By Joshua U. Otaigbe, PhD | March 17, 2026 | 6 min read

Flaney Associates | Materials Engineering & Innovation

Introduction: The Problem With Standard Sizes

Human anatomy does not come in standard sizes. Every person's bones are shaped differently, influenced by genetics, age, injury history, and disease. Yet for most of modern medicine's history, orthopedic implants have been manufactured in a limited range of standard sizes and shapes. Surgeons have had to select the closest available option and adapt the surgical procedure to make it fit, sometimes removing additional bone to accommodate a standardized implant.

Additive manufacturing is changing this equation fundamentally. By combining medical imaging data with advanced materials processing, surgeons and engineers can now design and print implants that match a specific patient's anatomy with sub-millimeter precision. This is not a futuristic concept. Patient-specific 3D-printed implants are already being used in thousands of surgeries annually, and the technology is advancing rapidly.

From Scan to Implant: The Workflow

The process of creating a patient-specific implant begins with high-resolution medical imaging, typically CT scanning that creates a three-dimensional map of the patient's anatomy. Engineers use specialized software to segment the images, identifying bone boundaries and defects, and design an implant that fills the defect precisely while accounting for the mechanical loads the implant will need to bear.

This digital design is then transferred to a metal additive manufacturing system that builds the implant layer by layer from titanium or cobalt-chrome powder. The entire process from CT scan to finished implant can be completed in days, compared to weeks or months for custom implants.

made through conventional manufacturing. For patients with complex tumor resections, traumatic bone loss, or unusual anatomy, this speed can be the difference between a straightforward recovery and an extensive, complicated surgical procedure.

Why Titanium?

Titanium alloys, particularly Ti-6Al-4V, have become the dominant material for orthopedic additive manufacturing, and for good reason. Titanium combines an excellent strength-to-weight ratio with exceptional corrosion resistance in the body's chemical environment. Critically, it is one of only a handful of metals that the human body accepts without a significant immune response, a property called osseointegration.

But additive manufacturing does not just replicate what could be done with conventionally manufactured titanium. It enables a design strategy that is impossible with any other manufacturing method: porous structures that mimic the architecture of natural bone.

"A solid titanium implant and a porous printed implant of the same size can behave completely differently inside the body. The architecture is as important as the material."

The Science of Osseointegration Through Porosity

Natural bone is not solid. It has a hierarchical porous structure, from the large cavities of trabecular bone visible to the naked eye, to microscopic channels and pores that allow blood vessels and bone-forming cells to permeate the entire structure. When bone grows into and around an implant, this process is called osseointegration, and it is what gives a permanent implant its long-term stability.

Smooth, solid metal implants rely primarily on surface chemistry and mechanical fixation for stability. Porous additive manufactured implants can provide a three-dimensional scaffold into which bone actively grows, creating a biological bond far stronger than any mechanical fixation alone. Research shows that pore sizes between 300 and 600 micrometers, with interconnected porosity of 65 to 80 percent, are optimal for bone ingrowth. Additive manufacturing can produce these structures with precise, reproducible geometry.

Regulatory Considerations for Printed Implants

The regulatory pathway for 3D-printed implants adds complexity to an already challenging field. The FDA classifies patient-specific 3D-printed implants as custom devices, which have a different

regulatory pathway than standard manufactured implants. While custom devices can be exempt from full premarket approval, the manufacturer must still demonstrate that the device is safe and effective for its intended use.

Critical regulatory requirements for additively manufactured implants include:

- Validation of the design software and the digital-to-physical manufacturing chain
- Demonstration that the material properties of printed parts meet the same standards as conventionally made implants
- Characterization of residual porosity, surface roughness, and subsurface defects
- Verification that post-processing steps (heat treatment, surface finishing, sterilization) do not degrade critical properties
- Traceability from patient imaging data through design, manufacturing, and implantation

The Future: Bioprinting and Smart Implants

The frontier of implant additive manufacturing extends beyond metals. Researchers are actively developing techniques for bioprinting, where living cells are printed into three-dimensional scaffolds to create tissues and ultimately organs. While fully bioprinted organs remain years away from clinical use, bioprinted cartilage patches, bone scaffolds, and vascular grafts are already in early clinical trials.

At Flaney Associates, we support medical device companies at the intersection of additive manufacturing, advanced materials, and regulatory compliance. Whether you are developing a patient-specific implant program or scaling an existing additive manufacturing process, our team provides the materials expertise to help you move confidently from concept to clinical use.

For more information or if you have any questions, please contact the author:

Joshua U. Otaigbe, PhD

Founder & Principal Consultant, Flaney Associates

Email: info@flaneyassociates.com | Web: flaneyassociates.com

[Schedule a free consultation at flaneyassociates.com/contact](https://flaneyassociates.com/contact)

© 2026 Flaney Associates. All rights reserved. This article is provided for informational purposes only.